

Adult Bibs

Disposable hygienic bibs for adults. The **Adult Bibs** is made of white absorbent paper and waterproof PE film.

It has a collection pocket at the bottom for good hygiene and ribbons to tie it at the back of the neck.



PRODUCT CHARACTERISTICS

- Colour: White paper and blue film
- Size 36x66 cm + 10 cm collection pocket
- Materials: paper layer 20 g/m² & film layer 15 g/m²



INSTRUCTIONS FOR USE

Turn the pocket around so that the white paper is on the outside.

Gently stretch the upper die-cut area to obtain two tapes.

Place the bib with the plastic area in contact with the user.

Discard after use.

CONFIGURATION

REFERENCE	DESCRIPTION	UNITIES BOX	MEASUREMENTS BOX (cm)	WEIGHT BOX* (kg)
8100000	ADULT BIBS 36x66+10 CM	500 Un. 4 bags of 125 Un.	38x36x15	4.8
8100000E	ADULT BIBS 36x66+10 CM	1000 Un. 8 bags of 125 Un.	40x40x20	10.7
8100000G	ADULT BIBS 36x66+10 CM	1000 Un. 20 bags of 50 Un.	40x40x20	10.7

* TOLERANCES: Weight ± 0.5 Kg