

# Eye Protector

## *Nursing Care – Paediatric Product*

### INTENDED USE

The **Eye Protector** is a medical device which is intended to be used as paediatric eye protector for phototherapy.

Easily adjustable thanks to the velcro system.

### COMPOSITION

Product with a three-layer material: cotton inner layer, outer polyamide layer and polyurethane foam filling.



### INSTRUCTIONS FOR USE

- Place the Protector around patient's eyes, around the head.
- Hold the Protector by means of the velcro system.

### PRODUCT CHARACTERISTICS

The product is adjustable by means of the velcro system, the elastic rubber strip and the material of the outer layer. The product is breathable, not irritating, antiallergic. The velcro fastening system allows an effective subjection.

### PRECAUTIONS AND RECOMMENDATIONS

The use by more than a patient may provoke a cross-infection.  
Do not use the product on injured skin.  
Do not use the product if it is polluted, dirty or evidently deteriorated.

### STORAGE AND SHELF-LIFE

Keep away from sunlight. Keep dry.  
Shelf-life: 5 years.

**PRODUCT DIMENSIONS**

In order to choose the correct size, the cranial perimeter has to be taken into account, as well as, the ocular surface to protect:

REFERENCE	MODEL	CRANIAL PERIMETER	OCULAR SURFACE	UPPER STRIP LENGTH
1029050	EXTRA-SMALL	25 – 30 cm	4 x 8 (cm)	25 cm
1029100	SMALL	30 – 35 cm	5 x 8 (cm)	29 cm
1029150	MEDIUM	32 – 37 cm	6 x 9 (cm)	35 cm

CONFIGURATION			
REFERENCE	UNITS/BOX	BOX DIMENSIONS (cm)	BOX WEIGHT (kg)
1029050	25	20x15x15	0.4
1029100	25	20x15x15	0.4
1029150	25	20x15x15	0.42

**BARCODE (box, EAN14):**

Reference: 1029050	18427835105196
Reference: 1029100	18427835101884
Reference: 1029150	18427835101891

**IN COMPLIANCE WITH:**

- Directive 93/42/EEC of the Council, Directive 2007/47/EC of the European Parliament and of the Council and Regulation (EU) 2017/745 of the European Parliament and of the Council.
- Class I, non-sterile (according to Annex IX of Directive 93/42 and Annex VIII of Regulation 745/2017).
- Enterprise certified by ISO 13485 and ISO 9001.
- Medical Devices Manufacturer Operating License No.: 3689-PS.

