

Operating Table Cover

INTENDED USE

The **Operating Table Cover** is a medical device, which purpose is to provide patients with the greatest of comfort during any operation or surgical procedure. Another use is to absorb the patient's blood or body fluids to clean the operating table or the patient better.

COMPOSITION

Bottom: Polypropylene 100% - waterproof.

Inside Filling: 100% Polyester Wadding synthetic high weight; anti allergic and antimould.

Upper: Non-woven fabric 100% viscose.

Latex free.



INSTRUCTIONS FOR USE

- Remove the product packaging and deploy the coating.
- Place the cover on the table. Make sure that the impermeable part is below.
- Discard after use.

PRODUCT CHARACTERISTICS

Helps maintain patient comfort during the course of the operation.

Bottom made of waterproof fabric to prevent going through the blood or body fluids to the surgical table. Single use.

PRECAUTIONS AND RECOMMENDATIONS

The use by more than one patient may cause a cross-infection.

Do not use the product on injured skin.

Do not use the product if it is pollutes, dirty or evidently deteriorated.

STORAGE AND SHELF-LIFE

Keep away from sunlight. Keep dry.

Shelf life: 5 years

MEDIDAS DE LA UNIDAD DE PRODUCTO

REFERENCE	MODEL	L (cm)	W (cm)
1040900	SMALL	66	23
1040910	MEDIUM	98	63
1040920	LARGE	194	69

CONFIGURATION			
REFERENCE	UNITS BOX	MEASUREMENTS BOX (cm)	WEIGHT BOX* (kg)
1040900	30 un/box	40x26x47	4.2
1040910	10 un/box	60x40x40	5.8
1040920	10 un/box	76.5x62x46.5	11.3

TOLERANCE: Weight $\pm 10\%$
IN COMPLIANCE WITH:

- Directive 93/42/EEC of the Council, Directive 2007/47/EC of the European Parliament and of the Council and Regulation (EU) 2017/745 of the European Parliament and of the Council.
- Class I, non-sterile (according to Annex IX of Directive 93/42 and Annex VIII of Regulation 745/2017).
- Enterprise certified by ISO 13485 and ISO 9001.
- Medical Devices Manufacturer Operating License No.: 3689-PS.



Sistema de Gestión
ISO 9001:2015
ISO 22716:2007
www.tuv.com
ID: S105678921