

Elbow Immobilizer

INTENDED USE

The **Elbow Immobilizer** is a medical device, which is intended to be use for the immobilization of the elbow, keep at arm's length and avoid bending.

It is placed directly on the arm, at the level of the elbow, a safe and comfortable result for the patient.



COMPOSITION

Inner part: Wood 100%.
Outer part: Polyester 100%.
Straps: Velcro

INSTRUCTIONS OF USE

- Place the immobilizer in the arm above the elbow.
- Secured by two side locks that incorporate Velcro system.

PRODUCT CHARACTERISTICS

Fast and easy to apply, adjustable and safe. Breathable, prevents skin irritation, suppresses overheating, moisture and maceration. Rounded edges.

PRECAUTIONS AND RECOMMENDATIONS

The use by more than one patient may cause a cross-infection.
Do not use the product on injured skin.
Do not use the product if it is pollutes, dirty or evidently deteriorated.

STORAGE AND SHELF-LIFE

Keep away from sunlight. Keep dry.
Shelf life: 5 years.

UNIT PRODUCT MEASUREMENTS

REFERENCE	MODEL	LENGHT (cm)	WIDTH (cm)
1060300	PEQUEÑO	14	14.5
1060310	GRANDE	19	19

CONFIGURATION			
REFERENCE	UNITS BOX	MEASUREMENTS BOX (cm)	WEIGHT BOX (kg)
1060300	24	21x15x15	0.536
1060310	24	21x15x15	0.890

TOLERANCE: Weight \pm 0.2 Kg

IN COMPLIANCE WITH:

- Directive 93/42/EEC of the Council, Directive 2007/47/EC of the European Parliament and of the Council and Regulation (EU) 2017/745 of the European Parliament and of the Council.
- Class I, non-sterile (according to Annex IX of Directive 93/42 and Annex VIII of Regulation 745/2017).
- Enterprise certified by ISO 13485 and ISO 9001.
- Medical Devices Manufacturer Operating License No.: 3689-PS.

