

Catheter holder

INDICATIONS

The **Catheter Holder** is indicated for a practical and efficient fixation of sonda-catheter to extremities.

COMPOSITION

Non-woven Polyester 100%
Velcro Polyamide 100%
Rubber band Polyamide 100%

INSTRUCTIONS FOR USE

- Close the velcro subsection.
- Place the tube and tie the velcro key

FEATURES

Designed with a large adjustable band and an elastic tissue to avoid the haemostatic effect. Made of a special extra soft tissue that combined with an adjustment system of Velcro provides a secure fixing and avoids the movement. Breathable material. Washable at 40°/60°C.

PRECAUTIONS AND RECOMMENDATIONS

The use by more than a patient may provoke a cross-infection.
Do not use the product on injured skin.
Do not use the product if it is polluted, dirty or evidently deteriorated.

STOREGE AND SHELF-LIFE

Keep away from sunlight. Keep dry.
Shelf life: 5 years.

PRODUCT DIMMENSIONS

REFERENCE	SIZE	L (cm)	W (cm)
1070570	STANDARD	62	5

PRESENTATION			
REFERENCE	UNITS / BOX	BOX DIMMENSIONS (cm)	BOX WEIGHT (kg)*
1070570	30	20x15x15	0.52

TOLERANCE: WEIGHT $\pm 10\%$

IN COMPLIANCE WITH:

- Directive 93/42/EEC of the Council, Directive 2007/47/EC of the European Parliament and of the Council and Regulation (EU) 2017/745 of the European Parliament and of the Council.
- Class I, non-sterile (according to Annex IX of Directive 93/42 and Annex VIII of Regulation 745/2017).
- Enterprise certified by ISO 13485 and ISO 9001.
- Medical Devices Manufacturer Operating License No.: 3689-PS.

