

# Monitor Set Holder

## INDICATIONS

**Monitor set holder** is indicated for the transport of equipment for monitoring purposes. Practical and comfortable. With internal protection. Standard dimensions.

## COMPOSITION

1070600-610: Non-Woven polyester 100%  
1070610: Polyethylene 100% - inner layers

## INSTRUCTIONS FOR USE

- Introduce the monitoring device into the pocket of the monitor set holder.
- Place the monitor set holder in the more useful and comfortable way, using the straps with adjustable velcro closure system.



## FEATURES

Manufactured with a special material with inside reinforcement for higher strength. The monitor set holder is perfectly adjusted to the body of patient with two long straps firmly fixed with a velcro closure system. Standard dimensions.

## PRECAUTIONS AND RECOMMENDATIONS

The use by more than a patient may provoke a cross-infection.  
Do not use the product on injured skin.  
Do not use the product if it is polluted, dirty or evidently deteriorated.

## STORAGE AND SELF-LIFE

Keep away from sunlight. Keep dry.  
Shelf life: 5 year

## PRODUCT DIMENSIONS

REFERENCE	LENGTH (cm)	WIDTH (cm)	TAPE LENGTH (cm)
1070600 1070610	20	16	120

PRESENTATION			
REFERENCE	UNITS/BOX	BOX DIMENSIONS (cm)	BOX WEIGHT (kg)*
1070600 1070610	25	21x15x15	0.71

TOLERANCE: ± 0.3 kg

## IN COMPLIANCE WITH:

- Directive 93/42/EEC of the Council, Directive 2007/47/EC of the European Parliament and of the Council and Regulation (EU) 2017/745 of the European Parliament and of the Council.
- Class I, non-sterile (according to Annex IX of Directive 93/42 and Annex VIII of Regulation 745/2017).
- Enterprise certified by ISO 13485 and ISO 9001.
- Medical Devices Manufacturer Operating License No.: 3689-PS.

