

# 3-Way Stopcock

*With extension*

## INTENDED USE

**3-way stopcock with extension** is a medical device sterilized using Ethylene Oxide gas and easy visualization to be able identified without opening.

Three ways, two female luer-lock connections and one male. With different length extension tubes and protective cap.

This product is suitable for lipids, cytostatics, antiseptics and other drugs of high caustic power.



## COMPOSITION

Polycarbonate body completely transparent. Extension tube: PVC DEHP FREE Ø 2.6x3.9 mm (8400022 Ø 3.0x4.1 mm). Exempt from latex and pyrogen.

## INSTRUCTIONS FOR USE

Recommended use for intravenously infusion and intermittent administration of medications.

- Remove the protective caps and connect the 3-way stopcock to the terminals.
- Make the liquid flow and move the air from the tube.
- Connect the connector of the extension cable with the corresponding terminal.

## PRODUCT CHARACTERISTICS

Product sterilized with Ethylene Oxide (ETO) gas as per standarized and validated sterilization cycle as per ISO 11135.

Compatible with blood and blood products.

Red and blue pegs for arterial and venous line identification.

Pressure tested pneumatically and hydrostatically: 6 bars.

## PRECAUTIONS AND RECOMMENDATIONS

Single use.

The use by more than one patient may cause a cross-infection.

Do not use the product if it is polluted, dirty or evidently deteriorated.

## STORAGE AND SHELF-LIFE

Keep away from sunlight. Keep dry.

Shelf-life: 5 years.

**PRODUCT DIMENSIONS**

REFERENCE	EXTENSION TUBE LENGTH (cm)
8400021	10
8400022	15
8400023	25
8400024	50
8400025	100

Packaged in unit blister with product description

CONFIGURATION			
REFERENCE	UNITS BOX	MESUREMENTS BOX (cm)	WEIGHT BOX* (kg)
8400021	200 Un. (8 bags of 25 un.)	39x38x29	2.49
8400022	200 Un. (10 bags of 20 un.)	39x39x28	2.58
8400023	200 Un. (10 bags of 20 un.)	39x39x28	2.69
8400024	200 Un. (10 bags of 20 un.)	39x39x28	3.10
8400025	200 Un. (10 bags of 20 un.)	39x39x28	3.97

\* TOLERANCE: Weight  $\pm$  0.5 Kg

**IN COMPLIANCE WITH:**

- Directive 93/42/EEC of the Council, Directive 2007/47/EC of the European Parliament and of the Council and Regulation (EU) 2017/745 of the European Parliament and of the Council.
- Class I, sterile (according to Annex IX of Directive 93/42 and Annex VIII of Regulation 2017/745).
- Enterprise certified by ISO 13485 and ISO 9001.
- Medical Devices Distributor License No.: 3689-PS.

