

# DISPOMEDIC DRY SCRUB

*Soft Scrub brush sterile for handwashing*

## DESCRIPTION

Disposable Dry sterile Brush/Sponge suitable for washing the hands of technical-healthcare personnel. Hand washing in high risk zones. Includes nail cleaner. Hygiene with preoperative patients.

## COMPOSITION

BRUSH: Polyethylene  
 SPONGE: 100% polyurethane foam  
 NAIL CLEANER: Copolymer ABS (acrylonitrile butadiene styrene)  
 Latex free



## INSTRUCTIONS

The purpose of this product is to clean the hands and nails of patients and healthcare personnel. It must be used as follows: Remove the brush from the wrapper, using the easy-open device located at one end of the pack. Once the brush has been removed from the wrapper, hold it by the middle and dispose of the wrapper. Clean the hands and nails, if necessary impregnating the sponge with an antiseptic, which is not included with this product. Use the nail cleaner for difficult access zones. Dispose of after use. Do not re-use.

## CHARACTERISTICS

Brush with bristles that holds a soft sponge. Its anatomical design ensures it is more convenient and flexible for the user. The nail cleaner is a thin green solid stick, easy to use and with high tensile strength and tear. Sterilised by ethylene oxide

## DIMENSIONS OF PRODUCT UNIT

80x50x37 mm. ±2.5 mm

PRESENTATION				
REFERENCE	DESCRIPTION	UNID/BOX	BOX MEASUREMENTS (cm)	BOX WEIGHT* (kg)
0030200F	DISPOMEDIC DRY SCRUB	100	38 x 21 x 35	1.8
0030200H	DISPOMEDIC DRY SCRUB 15 dispensers of 30 units	450	60 x 39.5 x 50.5 (dispenser 49.5x12.5x11.7)	10.4

## COMPLIES WITH:

-Quality management systems-requirements (ISO 9001:2015)

