

PRODUCT CATALOGUE

cvmédica®



Hygiene



Antisepsis



Nursing
care



Medical
devices



2023
2024

Index



HYGIENE

Pg. 5

Fibre sponges with soap	6
Foam sponges with soap	8
Foam soap-free sponges	9
Wash mitts with soap	10
Wet wash mitts	12
Wet wipes	13
Body hygiene liquids	14

ANTISEPSIS

Pg. 15

Surgical brushes with chlorhexidine	16
Antiseptic liquids with chlorhexidine	17
Antiseptic wash sponges	18
Brushes and dispensers	19
Waterless antiseptic wash mitts for skin	20



NURSING CARE

Pg. 21

Immobilization	22
Hand and paediatric immobilisers	23
Immobilisers and slings	24
Protection and anti-bedsore	27
Anti-bedsore cushions	28
Eye protection and coverings	29
Supports	30
Jockstraps	30

MEDICAL DEVICES AND PROTECTIVE PRODUCTS

Pg. 31

Gloves	32
Masks	33
Shields	34
Bibs	34
Syringes	35
3-part syringes	35
3-way stopcock and Infusion sets	36
Care for Incontinence and Urology	37
Disposable absorbent underpad	37
Medical Sheet	37
Continence Care and Urology	38
Non-sterile clothing	39
Single-use clothing	40
Accessories	42
Caps	42
PLUS caps	43
Patient kits	45
Disposan-tech reusable products	46
Reusable uniforms	48

ANNEX PRODUCTS, REFERENCES, LANGUAGES AND PACKAGING


Pg. 49

Please ask about product availability



Founded in 1991

Over **30 years** serving you.



Extensive portfolio


More than **500 references**.

cvmedica®

We are a specialist Spanish company engaged in the manufacture and distribution of healthcare products.

We provide solutions in hygiene, antisepsis, nursing care, protection equipment and medical devices intended for healthcare professionals, companies and users.

We began our journey with the manufacture of soapy sponges and a line of nursing products, growing until reaching more than 500 products, which has let us position ourselves in the hospital and geriatric field in more than 60 countries worldwide.



Manufacturing

More than **8,000m2** with highly automated manufacturing lines.



Exports

To more than **60 countries** since 1993.



Our team

Diverse, plural and trained, from more than **14 countries**.



Stock

More than **5,000m2** distributed over **two warehouses**.



Company policies



Satisfaction

Satisfy workers and the users of our products by improving their quality and safety.



Commitment

Consolidate and maintain relationships with stakeholders that are relevant to the growth of **CV Médica**.



Quality

Strictly comply with all the legal and regulatory requirements that apply not only to **CV Médica** as a company, but also to its products.




Constant improvement

Ensure that every employee is involved in improving our productivity, competitiveness, quality, prevention, asset management and respect for the environment.



Reliability

Maintain our customers' trust by satisfying their needs quickly and flexibly.



Equality

Our company is firmly committed to the well-being of our employees. As a result, it promotes measures such as equality and training plans that promote equal access to opportunities, and maintaining a suitable work environment.



Quality certificates



ISO 9001: 2015

Quality Management Systems.

ISO 13485: 2016

Medical devices. Quality Management Systems.

ISO 22716: 2008

Cosmetic Products. Good Manufacturing Practices (GMP).

Symbology



For general use



Suitable for babies



Nursing homes



Sports



Reduced mobility



Hospital



Travel



Industrial use



Dermoprotector



Dermatologically tested



Fresh aroma



No phthalates



No latex



PVC-free



With aloe vera extract




Biodegradable




100% cellulose

(*) Available in other formats. See annex with references.



Hygiene products

Full range of products for personal hygiene. The formulation of the products in this range respects the skin's neutral pH to help keep it in optimum condition.





Sponges

CV Médica manufactures two types of sponges - fibre and foam - both impregnated with dry soap with a neutral pH of 5.5.

Fibre sponges with soap



Dispobaño BIO

Dimensions: 20x12cm



Ref.	Units/Pack
200	10



Dispobaño Plus Aloe

Weight: 130g/m²

Dimensions: 20x12cm

With aloe vera extract.



Ref.	Units/Pack
112C	20



Dispobaño Aloe Vera

Weight: 100g/m²

Dimensions: 20x12 cm

With aloe vera extract.



Ref.	Units/Pack
104*	24

Instructions for use



Take a sponge and slightly dampen it.



Squeeze several times until foam appears.



Lather and rinse with water.



It is not recommended to immerse the sponge in water, as the soap can be lost.



Dispobaño

Weight: 100g/m²

Dimensions: 20x12 cm

Ref.	Units/Pack
100*	24



Dispobaño Plus

Weight: 150 g/m²

Dimensions: 20x12 cm

Ref.	Units/Pack
102*	20



Dispobaño Plus Compact

Weight: 170g/m²

Dimensions: 20x12 cm

Ref.	Units/Pack
106D	10

* See annex with references.



Foam sponges with soap



Dispofoam

Density: 15 Kg/m³

Dimensions: 20x12x1 cm

Ref.	Units/Pack
5403*	20



Dispofoam Aloe Vera

Density: 15 Kg/m³

Dimensions: 20x12x1 cm

With **aloe vera** extract to moisturise and soften skin.

Ref.	Units/Pack
5404*	20



Dispofoam One

Density: 15 Kg/m³

Dimensions: 13x8x2.5 cm

Ref.	Units/Pack
5412*	50



Dispofoam One Aloe Vera

Density: 15 Kg/m³

Dimensions: 13x8x2.5 cm

With **aloe vera** extract to moisturise and soften skin.

Ref.	Units/Pack
5405*	50

Instructions for use



Take a sponge and slightly dampen it.



Squeeze several times until foam appears.



Lather and rinse with water.



It is not recommended to immerse the sponge in water, as the soap can be lost.



Dispofoam One Plus

Density: 15 Kg/m³

Dimensions: 13x10x2.5 cm

Ref.	Units/Pack
5408	50

Foam soap-free sponges



Disposin One Plus

Density: 15 Kg/m³

Dimensions: 13x10x2.5 cm

Soap-free sponge.

Ref.	Units/Pack
5409	50



Disposin One

Density: 15 Kg/m³

Dimensions: 13x8x2.5 cm

Soap-free sponge.

Ref.	Units/Pack
5410*	50

* See annex with references.



Wash mitts with soap

- Anatomical shape
- Single 25x17cm size
- With non-woven and fibre fabric
- Impregnated with dry pH neutral (5.5) soap.
- Resistant and soft for the skin



Dispoglove Aloe Vera

Soap-impregnated fibre front and non-woven back side.

With aloe vera extract.

Ref.	Units/Pack	Units/Box
308	20	800



Dispoglove Protect

Soap-impregnated fibre front and non-woven back side with **waterproof interior.**

Ref.	Units/Pack	Units/Box
300	20	800



Dispoglove

Soap-impregnated fibre front and non-woven back side.

Ref.	Units/Pack	Units/Box
301	20	800

Instructions for use



Slightly dampen.



Lather and rinse with water.

It is not recommended to immerse in water, as the soap can be lost.



Disposin Glove Plus

- Fibre back side and front side
- Soap-free

Ref.	Units/Pack	Units/Box
303	20	600



Dispoglove Plus

- Soap-impregnated fibre front and fibre back side

Ref.	Units/Pack	Units/Box
307	20	600



Dispoglove Wash

- Front and back with non-woven fabric with waterproof interior, soap-impregnated front

Ref.	Units/Pack	Units/Box
309	50	2000



Instructions for use



Slightly dampen.



Lather and rinse with water.



It is not recommended to immerse in water, as the soap can be lost.



Disposin Glove Wash

- Front and back with non-woven fabric
- Soap-free

Ref.	Units/Pack	Units/Box
305C	50	2000



NEW

Disposin Mitt

- Front and back with non-woven fabric
- Soap-free

Ref.	Units/Pack	Units/Box
302	50	2500

Wet wash mitts

Wet mittens with non-soapy lotion.

- Made of non-woven fabric
- Full hygiene without water
- Ready to use
- Less time for personal hygiene
- No parabens, latex or fragrances
- Suitable for sensitive skin



Glowet

- Its composition allows it to be heated

Ref.	Units/Pack	Units/Box
408	8	288
412	12	288



Wet wipes

Wet wipes with non-soapy lotion

- Ready to use
- Safe for sensitive skin



Dispowipe Aloe

- Wet wipes with non-soap lotion
- Ready to use
- Enriched with Aloe Vera and chamomile extract

Ref.	Units/Pack	Units/Box
605	8	800
601	42	1050
602	80	1600

NEW

Dispowipe Aloe Bio

Ref.	Units/Pack	Units/Box
612	8	640
613	42	756
614	60	840



NEW

Dispowipe Baby

Ref.	Units/Pack	Units/Box
610	54	1296



NEW

Dispowipe Baby Bio

Ref.	Units/Pack	Units/Box
616	54	1296



NEW

Dispowipe Facial

Ref.	Units/Pack	Units/Box
617	10	800



NEW

Dispowipe Intimate

Ref.	Units/Pack	Units/Box
611	10	800



Body hygiene liquids

Cosmetics for full body hygiene.
Three products to clean, moisturise and scent the most delicate skin.



Dispobaño Gel-Shampoo 2 in 1

Bathroom and shower gel for daily hygiene.

Features:

- For all skin types
- Suitable for hair
- Contains moisturiser and moisturising agents
- No perfumes or dyes

Ref.	Size	Units/Box
700C	1000 ml with dispensing pump	12



Dispobaño Body Milk with Aloe Vera

Body milk that softens, nourishes and intensely moisturises the skin.

Features:

- Emollient and moisturising agents, such as aloe vera extract
- Suitable for sensitive and delicate skin

Ref.	Size	Units/Box
702C	1000 ml with dispensing pump	12



Dispobaño Eau de Cologne

Eau de cologne for daily personal hygiene.

Features:

- Low in alcohol
- Suitable for sensitive and delicate skin

Ref.	Size	Units/Box
701C	1000 ml	12





Antisepsis

Products to completely disinfect the hands and skin of healthcare professionals, patients and residents.



Chlorhexidine products

Surgical scrub brushes



Disposcrub C

Dimensions:

80x50x37 mm. + 2.5 mm.

Active ingredient:

With 4% chlorhexidine digluconate.

Ref.	Units/Box
30230*	100

Available with/without nail cleaner.



Dispomedic Scrub C Plus

Dimensions:

80x50x37 mm. + 2.5 mm.

Active ingredient:

With 0.8% chlorhexidine digluconate soap.

Ref.	Units/Box
30220*	100

Available with/without nail cleaner.

Antiseptic liquids

Products to properly disinfect the hands and skin of healthcare professionals, patients and residents.



Dispoclorhex 4 Soap Scrub

Use:

Desinfection of healthy skin.
Bactericide and yeasticidal/fungicide.

Active ingredient:

With 4% chlorhexidine digluconate soap.

Ref.	Content
40400*	500 ml* with dispensing pump



Dispomedic Soap 0,8%

Use:

Desinfection of healthy skin.
Bactericide and fungicide.

Active ingredient:

With 0.8% chlorhexidine digluconate soap.

Ref.	Content
40080*	500 ml* with dispensing pump



Dispoyod Povidona iodine 10%

Function:

Fungicide and bactericide /yeasticide.

Active ingredient:

10% Povidone-iodine solution.

Ref.	Content	Units/Box
DPY0050	50ml	96
DPY0125	125ml	12
DPY0500	500ml	12
DPY1000	1000ml	8



Dispomedic soap solution 7,5%

Function:

Fungicide and bactericide /yeasticide.


Active ingredient:

7.5% povidone iodine soap.

Ref.	Content	Units/Box
DPM0500	500ml with dosing pump	12
DPM1000	1000ml with dosing pump	8



* See annex with references.



Antiseptic wash sponges

Product for antiseptic skin cleaning.



Disposcrub Sponge C

Dimensions:

100x100x25 mm.

Active ingredient:

20 ml of 4% chlorhexidine digliconate soap.

Ref.	Units/Box
30330*	120



Dispomedic Sponge C Plus

Dimensions:

100x100x25 mm.

Active ingredient:

20 ml of soap with 0.8% chlorhexidine digluconate soap.

Ref.	Units/Box
30320	120

Brushes and dispensers

Brush for cleaning the hands and nails of technical-sanitary personnel.



Dispomedic Scrub Dry

- Sterilised with ethylene oxide

Dimensions:

80x50x37 mm.

Ref.	Units/Box
30200*	100



Dispomedic P Scrub

Dimensions:

80x50x37 mm.

Active ingredient:

20ml of 7.5% iodopovidone soap.

Ref.	Units/Box
30210*	100



Metallic dispensers

- Metal dispensers for brushes or brush dispensers
- Made of AISI 304 stainless steel
- Easy to clean and suitable for operating theatres

Ref.	Stand for	Dimensions	Units/Box
30007	Dispensing box 30 units	13x13x16 cm.	1
30009	Rushes 15 units	17.4x13x60 cm.	1

* See annex with references.

Waterless antiseptic wash mitts for skin



Glowet dispoclorh

Mitten with 2% colourless aqueous chlorhexidine.

Use:

- Helps to reduce bacteria that cause skin infections
- Recommended for pre-operative bathing and complete hygiene of patients with bedsores or reduced mobility
- One-size-fits-all: 25x17 cm
- Anatomical shape

Ref.	Units/Pack	Units/Box
40010	6	240
40011	8	288



Nursing care

Products that offer the necessary functionality, with unbeatable efficiency for healthcare professionals and maximum comfort for the user.





Immobilization

Limb immobiliser

For restless patients who risk stopping treatment, injuring themselves or jeopardising the integrity of the affected area. With an anatomical design tailored to the wrists and ankles, it includes a closing system to adjust the immobiliser to the extremity with strips that can be attached to the bed.

Limb holder with Orifice

Material: Non-woven fabric with foam interior.

Colour: White.

Features:

- Strong fastening through hole

Ref.	Units/Box
1090303	80



Fast Opening Limb Holder

Material: Non-woven fabric with foam interior.

Features:

- With velcro closure for quick opening

Ref.	Units/Box
1090110	60



Safety Limb Holder

Material: Non-woven fabric with foam interior.

Features:

- With velcro closure and ring for greater safety

Ref.	Units/Box
1090100	80



Comfort Safety Limb Holder

Material: Cotton interior and padded exterior.


Features:

- Velcro closure and ring for greater safety

Ref.	Sizes	Units/Box
1000202	L	60
1000405	S	60



Hand and paediatric immobilisers



Limb Holder with Hand Protection Glove

Material: Breathable mesh on the back of the hand and soft, padded fabric for the palm.

Features:

- Recommended for hand surgery and immobilisation of patients
- Provides finger clearance and offering the best posture for the hand and its immobilisation

Ref.	Units/Box
1090300	16



Pediatric-Neonatal Limb Holder

- Padded material, very soft on the skin, with velcro fastener

Ref.	Units/Box
1090150	120

Elbow immobiliser for children

Material: Exterior in non-woven fabric, interior with wooden depressor.

Features:

- Comfortable, lightweight and breathable
- Adjustable with velcro closures



Ref.	Sizes	Units/Box
1060310	M	24
1060300	S	24



Immobilisers and slings

Arm and shoulder immobiliser: restrains the arm and shoulder with an adjustable thoracic strip.

The large size of the pouch for the forearm holds the wrist, includes a thumb strap for maximum comfort. It can be used as a sling.

Sling: supports the patient's arm, whether is fractured or with minor injuries in recovery phase.



Grey shoulder immobiliser

Material: Padded material with finger support and adjustable straps for added comfort.

Ref.	Sizes	Units/Box
1000720	XL	10
1000721	L	10
1000722	M	10
1000723	S	10
1000724	XS	10



Blue shoulder immobiliser

Material: Lightweight, breathable, non-woven fabric with padded straps for added comfort.

Ref.	Sizes	Units/Box
1012915	XL	20
1012925	L	20
1012935	M	20
1012945	S	20
1012955	XS	20



Mesh shoulder immobiliser

Material: Breathable material and adjustable straps in non-woven fabric.

Ref.	Sizes	Units/Box
1012825	L	10
1012835	M	10
1012845	S	10



Neck/cuff sling

Material: Made with a lightweight, breathable material. Velcro adjustment system.

Ref.	Sizes	Units/Box
1000710	Adulto	30
1000709	Pediátrico	30



Clavicle immobiliser

Holds the shoulder in the proper position to recover from a collarbone injury.

Material: Padded armpit straps and 4-point adjustment for greater comfort.

Ref.	Sizes	Units/Box
1000705	L	10
1000706	M	10
1000707	S	10
1000711	XS	10



Abdominal binder

Post-operative abdominal support.

Ref.	Sizes	Units/Box
1070310	L	10
1070300	M	10



Moldable support

Designed to immobilise and protect puncture sites on patients.

Ref.	Sizes	Units/Box
1060130	L	30
1060125	M	30
1060120	S	30
1060100	XS	30





IMPROVED PRODUCT

Comfort perineal fastening belt

Breathable perineal belt for perfect chair support.

Ref.	Units/Box
1070210	12
1070200	12



Normal support vest

Lightweight, highly resistant breathable material.

Features:

- Efficient and secure support

Ref.	Sizes	Units/Box
1015670	L	10
1015660	M	10
1015650	S	10



Abdominal support belt

Foam abdominal belt for bed or chair support.

Ref.	Sizes	Units/Box
1070150	L	12
1070100	M	12



Protection and anti-bedsore

Protective elements to prevent the appearance of bedsores in patients who have to undergo prolonged periods of immobilisation or reduced mobility.



Comfort elbow-heel protector

Features:

- Maximum comfort due to its foam padding and cotton lining

Ref.	Units/Box
1000610	60



Heel-foot protector

Features:

- Full protection of the foot, protects the heel, with non-woven fabric and fibre material
- Breathable

Ref.	Units/Box
1080250	20



Elbow heel protector

Features:

- Its polyester fibre filling distributes pressure evenly, protecting the articulation

Ref.	Characteristics	Units/Box
1080200	Standard	36
1080210	Reinforced	36



Open elbow-heel protector

Features:

- The low-friction resistance coating makes it difficult to remove the protector

Ref.	Dimensions	Units/Box
1080290	Short	40
1080280	Long	20

Anti-bedsore cushions



Dotted anti-bedsore cushion

The dotted anti-bedsore cushion is used to prevent bedsores in patients who have to stay in a chair for a long time.

Features:

- Made of cotton and silicone fibre. Comfortable and breathable

Ref.	Dimensions	Units/Bag
5020100	48x48 cm	1



Anti-bedsore

Cushion covered in synthetic soft fibre to prevent bedsores in patients who have to remain seated for a long time.

Available in different shapes that adapt to the needs of the patient or chair type.

Features:

- Made of cotton and siliconised fibre. Comfortable and breathable

Ref.	Characteristics/ dimensions	Units/Bag
5010001	Horseshoe (44x44 cm)	1
5010003	Ring (Ø 49 cm)	1

Eye protection and coverings



Eye protector

Eye protection for phototherapy to prevent possible injury in paediatric and newborn patients.

Features:

- Adjustable to cranial perimeter and head height

Ref.	Sizes/cranial perimeter	Units/Box
1029150	M 32-37 cm	25
1029100	S 30-35 cm	25
1029050	XS 25-30 cm	25



Eye protector with rear attachment

Eye protection for phototherapy to prevent possible injury in paediatric and newborn patients, with rear fastener.

Features:

- Adjustable to cranial perimeter and head height

Ref.	Sizes/cranial perimeter	Units/Box
1029220	L 35-41 cm	25
1029200	M 36-40 cm	25
1029250	S 28-35 cm	25



Braun splint cover

Comfortable, breathable lining to pad the Braun splint cover.

Features:

- Fibre interior

Ref.	Sizes	Units/Box
1040350	XL 122x68 cm	12
1040310	L 105x54 cm	12
1040305	M 95x54 cm	12
1000321	M 95x54 cm	12
1040300	S 85x54 cm	12



Surgical table overlay

Features:

- Breathable, comfortable and fluid-absorbent

Ref.	Sizes	Units/Box
1040920	L	10
1040910	M	10
1040900	S	30

Supports

Fastening elements to meet needs in an efficient and functional manner.



Catheter-probe support

This support is recommended for securely fastening probes and catheters to the limbs.

Features:

- Practical and easy to use

Ref.	Units/Box
1070570	30



Telemetry monitor holder

Suitable for transporting telemetry devices and monitors. Practical and convenient system.

Features:

- Standard measurements

Ref.	Units/Box
1070600	25

Jockstraps



Medical jockstrap

- Therapeutic use
- Soft and fully adjustable

Ref.	Units/Box
1070811*	10



Universal jockstrap

Jockstraps for the male genital area.

Features:

- Made with highly resistant elastic material
- Effective for therapeutic and prophylactic support
- Available in several sizes

Ref.	Units/Box
1070803*	10



Medical devices and protective products

Single-use disposable products.





Gloves

Disposable examination gloves certified as personal protective equipment (PPE).

- Suitable for handling food
- Powder free
- AQL 1.5



Nitrile gloves

Ref.	Size	Pcs./Pack	Units/Box
9000000	XS	100	1000
9000001	S	100	1000
9000002	M	100	1000
9000003	L	100	1000
9000004	XL	100	1000



Latex gloves

Ref.	Size	Pcs./Pack	Units/Box
9000010	XS	100	1000
9000011	S	100	1000
9000012	M	100	1000
9000013	L	100	1000
9000014	XL	100	1000



Vinyl gloves

Ref.	Size	Pcs./Pack	Units/Box
9000020	XS	100	1000
9000021	S	100	1000
9000022	M	100	1000
9000023	L	100	1000
9000024	XL	100	1000



Masks

Disposable surgical mask with elastic rubber or straps, 3 layers, blue.

Instructions for use



NEW

Dispomedic mask with straps

- Type IIR
- Mouldable nasal clip for better fit
- Fluid resistant
- Latex free
- Bacterial filtration efficiency > 98%

Ref.	Units/Pack	Units/Box
9000101	50	1000



Dispomedic mask with elastic band

- Type IIR
- Mouldable nasal clip for better fit
- Fluid resistant
- Latex free
- Bacterial filtration efficiency > 98%

Ref.	Units/Pack	Units/Box
9000100	50	2000

Mascarillas FFP2



Dispomedic FFP2 mask

- Respiratory protection mask
- Latex-free
- Integrated adjustable nasal clip
- Individually wrapped

Ref.	Units/Pack	Units/Box
9000200	50	1500

Shields



NEW

Face Shield with anti-fog

- To prevent contamination caused by liquid splashes or aerosols
- The headband has a foam base that relieves pressure
- It features an anti-fog and anti-static treatment

Dimensions: 33x22 ± 0,5

Ref.	Units/Pack	Units/Box
8880001	10	200

Bibs



Disposable bibs

Hygienic and disposable bib for adults. The adult bib is made of absorbent white paper and an impermeable PE film. With a collection pocket at the bottom for better hygiene, and straps to tie it to the back of the neck.

Dimensions

36x66 + 10 cm pocket

Ref.	Units/Pack	Units/Box
8100000*	125	500



Reusable bibs

Reusable adult bib with two layers. Made with an absorbent TNT layer and PVC backing, resistant to the passage of liquids.

Note:

- Size: 46x92 cm
- Absorption capacity: 800 ml
- Can be washed more than 150 times.
Do not use bleach

Ref.	Units/Box
8100100	50



Syringes

3-part syringes

- Made of non-toxic and apyrogenic polypropylene, and sterilised with ethylene oxide
- Easy to use



Luer cone syringes

Ref.	Tip	ml	Units/Pack	Units/Box
8500001	Luer	1 ml	100	3200
8500003	Luer	2.5 ml	100	3000
8500005	Luer	5 ml	100	2400
8500010	Luer	10 ml	100	1600
8500020	Luer	20 ml	50	800
8500050	Luer	50 ml	25	400



Luer lock cone syringes

Ref.	Tip	ml	Units/Pack	Units/Box
8501002	Luer Lock	2.5 ml	100	3000
8501005	Luer Lock	5 ml	100	2400
8501010	Luer Lock	10 ml	100	1600
8501020	Luer Lock	20 ml	50	800
8501050	Luer Lock	50 ml	25	400



Catheter cone syringes

Ref.	Tip	ml	Units/Pack	Units/Box
8501051	Catheter	50 ml	25	400
8500100	Catheter	100 ml	25	100

3-way stopcock and Infusion sets

Non-toxic and non-pyrogenic infusion sets after ethylene oxide sterilisation.

Recommended for intermittent infusion and the administration of intravenous medications.



Three-way stopcock

Three-way stopcock with two luer female and one male luer lock for controlling two ports of infusion/drainage through a single lancing device or line.

Features:

- Resistant to hydrostatic pressure
- Fully transparent polycarbonate body
- The product is suitable for lipids

Ref.	Units/Pack	Units/Box
8400020	50	500



Infusion therapy set

Equipment for parenteral administration of solutions or medications.

Ref.	Type	Units/Box
8400030	With air filter	200
8400031	With air filter and Y-shape injection point	200



Three-way stopcock with luer lock extension

Extender with three-way valve and two luer female and one male closures to control two infusion/drainage ports using a single lancing device.

Features:

- Resistant to hydrostatic pressure
- The product is suitable for lipids

Ref.	Dim.	Units/Pack	Units/Box
8400021	10 cm	20	200
8400023	25 cm	20	200
8400024	50 cm	20	200
8400025	100 cm	20	200



Three-way stopcock with rotating luer lock extension

Ref.	Extension length	Units/Pack	Units/Box
8400022	10 cm	25	200
8400026	25 cm	20	200
8400027	50 cm	20	200
8400028	100 cm	20	200



Care for Incontinence and Urology



NEW

Disposable absorbent underpad

- Prevents irritation or other lesions at skin level in patients
- High absorption capacity
- Distribution of fluids throughout the entire surface
- Breathable
- Latex and PVC free
- With a pattern that distributes fluids throughout the surface



Dispomedic underpad

Ref.	Dim.	Units/Pack	Units/Box
8880440	60x90cm	20	200
8880441	40x60cm	20	200
8880443	60x60cm	20	200

Dispomedic underpad PLUS

- Heavier

Ref.	Dim.	Units/Pack	Units/Box
8880442	60x90cm	20	160

Medical Sheet



Medical sheet

Surgical table pad.

Ref.	Dim.	Units/Pack	Units/Box
8880405	102x228cm	1	40
8880407	105x155cm	5	40

Medical sheet combo

Surgical table pad with patient cover.

Ref.	Dim.	Units/Pack	Units/Box
8880406	102x228cm	5	40



Contenance Care and Urology



Urine meter bag

- For draining urine during surgery, collecting urine or accurately measuring the urine volume
- Features a check valve and sampling port without needle

Measurement accuracy: ± 1 ml

Ref.	Units/Pack	Units/Box
7100010	1	10



Urine meter exchange bags

- For storing the urometer's urine when it requires emptying
- It features a holder for the drain tap

Ref.	Type	Units/Pack	Units/Box
7100011	Non-sterile	120	120
7100012	Sterile	1	60

Non-sterile cover and accessories





Single-use clothing

Non-sterile clothing for sanitary and industrial environments that require good manufacturing practices.

- Disposable, single-use products
- For full body coverage
- Comfortable and breathable



Isolation gown, green

- Single size
- Adjustable elasticated cuffs
- Adjustment with straps on the neck and waist

Ref.	Units/Pack	Units/Box
8200410	10	200



Isolation gown plus, green

- Single size
- Adjustable elasticated cuffs
- Adjustment with straps on the neck and waist

Ref.	Units/Pack	Units/Box
8201410	10	200



Isolation gown with fabric cuff, blue

- Single size
- 5-cm adjustable elasticated cuffs
- Adjustment with straps on the neck and waist

Ref.	Units/Pack	Units/Box
8200450	10	100



Isolation waterproof gown with fabric cuff, blue

- Waterproof
- Single size
- 5-cm adjustable elasticated cuffs
- Adjustment with straps on the neck and waist

Ref.	Units/Pack	Units/Box
8200402	10	100



Patient gown with short sleeves plus, blue

- Single size
- Non-transparent
- Short sleeves
- Adjustment with straps on the neck and waist

Ref.	Units/Pack	Units/Box
8200470	10	100



Patient gown with short sleeves plus, blue

- Single size
- Non-transparent
- Sleeveless
- Adjustment with straps on the neck and waist

Ref.	Units/Pack	Units/Box
8200460	10	100



Isolation gown, with Velcro front closure

- Single size
- Breathable
- 5-cm adjustable elasticated cuffs
- Velcro front closure

Ref.	Units/Pack	Units/Box
8200400	10	100



White protective Overall

- Zip fastening
- With hood
- Adjustable elasticated cuffs and hood

Ref.	Size	Dimensions	Units/Pack	Units/Box
8200500	S	158x54±2	1	50
8200501	M	162x56±2	1	50
8200502	L	162x56±2	1	50





Accessories

Complementary clothing products that cover the limbs and specific sections of the body.

Caps

- A garment that covers the hair completely and prevents it from falling into the environment
- SINGLE-USE product
- Odourless, breathable, hypoallergenic, anti-static and fire-retardant



Woman's bouffant cap, beret

- Single size
- Breathable

Ref.	Colour	Units/Pack	Units/Box
8200100	Blue	100	1000
8200110	Green	100	1000



NEW

Round nurse cap, green

Ref.	Units/Pack	Units/Box
8880010	100	2000



NEW

Surgeon cap with ties, green

Ref.	Units/Pack	Units/Box
8880015	100	2000

NUEVO

PLUS caps

- With dispensing box
- Equipment to fully cover the hair and prevent it from falling in the environment
- SINGLE-USE product
- Odourless, breathable, hypoallergenic, anti-static and fire-retardant



Special (bouffant) cap with overlock seam

- With elastic around the entire perimeter
- Breathable

Ref.	Units/Pack	Units/Box
8881010	250	1000



Surgeon cap plus with ties, green

- With straps at the back for improved fastening
- Breathable

Ref.	Units/Pack	Units/Box
8881015	150	600



Surgical cap with elastic

- Breathable
- With rubber at the back to hold long hair

Ref.	Units/Pack	Units/Box
8881018	250	1000



Cap astronaut type

- Full-head hood
- Covers the head and neck
- Breathable

Ref.	Units/Pack	Units/Box
8881021	250	1000





Non-woven fabric shoe cover

- Single size
- Elastic trim for better fit
- Breathable

Ref.	Colour	Units/Pack	Units/Box
8200210	Green	100	2000
8200200	Blue	100	2000



Shoe covers, green

- Single size
- Waterproof
- Elastic for better fit

Ref.	Units/Pack	Units/Box
8204210	100	2000



Blue sleeve covers

- With elastic for better fit

Ref.	Units/Pack	Units/Box
8200300	100	2000



Non-woven fabric beard snood, white

- Breathable
- Rubber trim

Ref.	Units/Pack	Units/Box
8200600	100	1000



Patient kits

Patient kit Private

Ref 8880203

- Isolation gown with cuffs
- Woman's cap, green, beret PLUS
- Panties/Slip
- Disposable slippers
- Toilet bag



Patient kit Outpatient surgery

Ref 8880201

- Patient gown braces PLUS
- Woman's cap, green, beret PLUS
- Shoe covers
- Panties/Slip



Patient Kit Consultations

Ref 8880202

- Patient gown braces PLUS
- Woman's cap, green, beret PLUS
- Shoe covers





**Disposan-tech
reusable
products**

Products developed with smart, reusable fabrics, designed mainly for the healthcare and geriatric sectors.



Sanitary Blanket

- Special blanket for operating theatre areas
- Suitable for home and healthcare use
- 3 sizes

Ref.	Size	Dimensions	Units/Box
8880402	S	72x80 cm	5
8880401	M	100x145 cm	5
8880400	L	145x200 cm	5



Reusable bib plus

- The fastening method prevents pulling from the front and makes it easy to open and close
- Adjustable to two measurements

Ref.	Units/Box
8100101	20



Draw sheet - Reusable underpad

- One piece, prevents breakage
- Uniform absorption across the entire surface

Ref.	Units/Box
8880430	5



NEW

Reusable uniforms

- Highly breathable garments, do not retain moisture
- Maintain the body's normal temperature
- Environmentally friendly, 100% recyclable



Warm-up scrub jacket plus

- Bomber jacket
- Fastening with pressure buttons
- With two front pockets

Ref.	Size	Units/Box
8880700	XS	5
8880701	S	5
8880702	M	5
8880703	L	5
8880704	XL	5
8880705	XXL	5



Basic cold scrub jacket

- Bomber jacket
- Fastening with pressure buttons
- With two front pockets

Ref.	Size	Units/Box
8880710	XS	5
8880711	S	5
8880712	M	5
8880713	L	5
8880714	XL	5
8880715	XXL	5



Scrub jacket

- Closed garment
- V-neck
- With three front pockets

COLOUR



Size	Ref.	Ref.	Ref.	Ref.	Units/Box
XS	8880500	8880510	8880520	8880530	5
S	8880501	8880511	8880521	8880531	5
M	8880502	8880512	8880522	8880532	5
L	8880503	8880513	8880523	8880533	5
XL	8880504	8880514	8880524	8880534	5
2XL	8880505	8880515	8880525	8880535	5
3XL	8880506	8880516	8880526	8880536	5



Scrub pants

- Unisex, adaptable to the body
- With adjustable drawstring
- Side pocket

COLOUR::



Size	Ref.	Ref.	Ref.	Ref.	Units/Box
XS	8880600	8880610	8880620	8880630	5
S	8880601	8880611	8880621	8880631	5
M	8880602	8880612	8880622	8880632	5
L	8880603	8880613	8880623	8880633	5
XL	8880604	8880614	8880624	8880634	5
2XL	8880605	8880615	8880625	8880635	5
3XL	8880606	8880616	8880626	8880636	5

ANNEX Products, references, languages and packaging

Hygiene

Dispobaño Aloe Vera

Ref.	Units/Pack	Units/Box	Tablet
104	24	3600	YES
2020104	24	1680	YES
2020104F	24	960	NO

Dispobaño

Ref.	Units/Pack	Units/Box	Tablet
100R	24	3600	YES
100S	24	1800	NO
100L	24	2400	YES
2020100R	24	1680	YES
2020100S	24	960	NO

Dispobaño Plus

Ref.	Units/Pack	Units/Box	Tablet
102	20	1200	NO
102B	20	2400	YES
102E	20	2000	YES

Dispofoam

Ref.	Units/Pack	Units/Box	Tablet
5403	20	2400	YES
2025403	20	800	YES

Dispofoam Aloe Vera

Ref.	Units/Pack	Units/Box	Tablet
5404A	20	2400	YES
2025404A	20	800	YES

Dispofoam One

Ref.	Units/Pack	Units/Box	Tablet
5412B	50	2000	YES
5412D	50	500	NO
5412E	50	1000	YES
5412H (MP)	50	2000	YES

Dispofoam One Aloe Vera

Ref.	Units/Pack	Units/Box	Tablet
5405A	50	2000	YES
5405B	50	500	NO

Disposin One

Ref.	Units/Pack	Units/Box	Tablet
5410D	50	2000	YES
5410E	10	500	NO
5410F (MP)	50	2000	YES

Antiseptis

Disposcrub C

Ref.	Nail cleaner	30-unit dispenser	Units/Box	Languages
30230	Without	NO	100	Spanish
30230B	With	NO	100	Spanish
30230C	Without	YES	450	Spanish
30230D	With	YES	450	Spanish
30230F	With	NO	100	English
30230G	Without	YES	450	English
30230H	With	YES	450	English

Dispomedic Scrub C Plus

Ref.	Nail cleaner	30-unit dispenser	Units/Box	Languages
30220	Without	NO	100	Spanish
30220B	With	NO	100	Spanish
30220C	Without	YES	450	Spanish
30220D	With	YES	450	Spanish

Dispoclorhex 4 Soap Scrub

Ref.	Content	Units/Box	Languages
40400	500 ml with dispensing pump	20	English
40400B	500 ml with dispensing pump	20	English
40400C	500 ml with dispensing pump	20	Portuguese
40400D	500 ml with dispensing pump	20	French
40410	1000 ml with dispensing pump	10	Spanish
40410B	1000 ml with dispensing pump	10	English
40410C	1000 ml with dispensing pump	10	Portuguese
40410D	1000 ml with dispensing pump	10	French

Dispomedic Soap 0.8%

Ref.	Content	Units/Box	Languages
40080	500 ml with dispensing pump	20	Spanish
40080B	500 ml with dispensing pump	20	English

Dispoclорhex 4 Soap Scrub (single-dose envelopes)*

* Minimum order on request.

Ref.	Content	Units/Box	Languages
40420	20 ml envelope	150	Spanish
40420B	20 ml envelope	150	English

DispoPET

Ref.	Content	Units/Box	Languages
40500	500 ml blind plug	20	Spanish
40500B	500 ml blind plug	20	English

Disposcrub Sponge C

Ref.	Units/Pack	Units/Box	Languages
30330	1	120	Spanish
30330B	1	120	English
30330C	1	120	Portuguese

Dispomedic Scrub Dry

Ref.	Nail cleaner	30-unit dispenser	Units/Box	Languages
30200	Without	NO	100	Spanish
30200B	With	NO	100	Spanish
30200C	Without	YES	450	Spanish
30200D	With	YES	450	Spanish
30200F	With	NO	100	English
30200H	With	YES	450	English

Dispomedic P Scrub

Ref.	Nail cleaner	30-unit dispenser	Units/Box	Languages
30210	Without	NO	100	Spanish
30210B	With	NO	100	Spanish
30210C	Without	YES	450	Spanish
30210D	With	YES	450	Spanish
30210F	With	NO	100	English
30210G	Without	YES	450	English
30210H	With	YES	450	English

Nursing care

Medical jockstrap

Ref.	Size	Dimensions	Units/Box
1070814	T4	12.75x11.5 cm	10
1070815	T5	13.75x12 cm	10
1070816	T6	14.5x13 cm	10
1070817	T7	15x13.5 cm	10
1070818	T8	15.75x14.5 cm	10
1070819	T9	16.5x15 cm	10
1070820	T10	17.75x15.75 cm	10
1070821	T11	18.5x16.75 cm	10
1070822	T12	19.5x17.25 cm	10

Universal jockstrap

Ref.	Size	Dimensions	Units/Box
1070803	T3	71/78 cm	10
1070804	T4	79/86 cm	10
1070805	T5	87/96 cm	10
1070806	T6	71/78 cm	10

Medical device

Disposable bibs

Ref.	Units/Pack	Units/Box
8100000	125	500
8100000E	125	1000
8100000G	50	1000

Ask about product availability



C.V. Médica, S.L. • C/Ponent, 1 – Pol. Ind. La Plana • 43424 Sarral, Tarragona SPAIN
Tel. +34 977 890 598 • info@cvmedica.com

www.cvmedica.com

Business card

cvmédica®

Expert manufacturers of
healthcare solutions



C.V. Médica, S.L. • C/Ponent, 1 – Pol. Ind. La Plana • 43424 Sarral, Tarragona SPAIN
Tel. +34 977 890 598 • info@cvmedica.com

www.cvmedica.com

